

Why Validation has become an Issue

DIY

Do-it-yourself crystallography is now a reality.

Instruments like the Bruker-Nonius kCCD and the Collect-Denzo-Scalepack suite of programs make reliable X-ray data collection quite practical with very little training.

Why Validation has become an Issue

Structure solution and refinement programs have also become less difficult to use.

In fact, SIR2002 and SHELXS are now so refined, there is probably very little the casual user can do if they fail to solve a structure.

The chemist now has a real possibility of determining their own crystal structures.

Why Do-it-Yourself Analyses?

Economics

We estimate that, excluding floor charges, depreciation etc, the actual running cost of a kCCD diffractometer is 12,000 Euro p.a.

This covers:

Electricity, Liquid Nitrogen, X-ray tube, routine maintenance, air conditioning, chilled water.

Why Do-it-Yourself Analyses?

Economics

In Oxford, except for special studies, the machines are run on a 3-shift system:

9.00-13.00 4 hours

13.00-1700 4 hours

Overnight 16 hours

This gives a capacity of over 900 crystal examinations pa, i.e. 13 Euro per experiment.

Why Do-it-Yourself Analyses?

Economics

In the UK, a young research assistant earns about 35,000 Euro pa.

Given a mix of easy and difficult organic & organo-metallic structures, a sustainable throughput might be 200 experiments pa, making the man-power cost per structure 175 Euro.

Why Do-it-Yourself Analyses?

Economics

Thus we have an instrument with a capacity of 900 investigations pa, and a research assistant with a capacity of 200 investigations pa. What is to be done with the remaining time?

A common solution is to reduce the throughput to only one investigation each day.

DIY Crystallography – Good or Evil?

The Good

- Improve chemists awareness of power and limitations of X-ray analysis.
- Make best use of crystallographers time.
- Make efficient use of diffractometers.

The Bad

- Inexperienced analysts can make mistakes.
- Careless users could wreck machine.

Chemical Validation

The aim of most structure analyses is to investigate physical or chemical problems. At the end of an analysis the questions to be asked are:

Is there any significant possibility that the structure may be 'wrong'?

What confidence can we have in the molecular parameters?

What confidence can we have in the adps?

Chemical Validation

If a structure ‘looks wrong’, it probably *is* wrong.

``X-ray analysis is still an art, through which our way
is lighted by crystal-chemical knowledge, and
chemical intuition.'`

The Cambridge Crystallographic Data Base stores much of our knowledge of what ‘looks right’, and MOGUL is a tool for accessing it.



What is MOGUL?

A data base of molecular geometries.


For every distinguishable chemical fragment in the CDS, the mean, standard deviation and distribution is computed.

Because the CSD has been pre-processed, MOGUL has super fast search speeds (a few seconds) by avoiding graph matching.



How does it do that?

By building a library of chemical fragments.






How does it do that?

A hierarchical key is built up

- element symbol
- number of bonded atoms
- hydrogen count
- *e.g.* C.3.1

- bond types
- *e.g.* 2.1


- Ring size *e.g.* 5





How does it do that?


- Search tree created



- Each leaf of the tree is searched for in CSD
- When matches are found, statistical quantities are computed and links set to actual structures.

How does it do that?

Keys are assigned to the user's query in the same way, and used to search the tree




If a match is found, statistical quantities are output and links set to actual structures.



MOGUL and CRYSTALS

Simple interaction. MOGUL Geometry
Check validates the whole structure

Richard Cooper
www.xtl.ox.ac.uk






MOGUL and CRYSTALS

Automatic interaction:

- Find all unusual bonds and valence angles.





MOGUL and CRYSTALS

Simple interaction – just click on the menu
MOGUL Geometry Check validates the
whole structure

Automatic interaction - find unusual bonds
and valence angles.

Each bond length and angle is assigned a z -
score:

$$z = \frac{\textit{reported} - \textit{CSDideal}}{\sigma(\textit{CSDideal})}$$



MOGUL and CRYSTALS

Automatic interaction - find unusual bonds and valence angles.

Mogul Results

Results | Atoms | Z Distribution | Information

Z-scores: zero is best agreement. See information tab.
Click a column header to sort data.

Query	Delta	Σ	Type
C(8) C(9) C(7)	-3.617	1.009	Angle
C(16) C(17) C(15)	-3.419	0.810	Angle
C(4) C(5) C(3)	-2.681	0.717	Angle
C(4) C(5)	0.030	0.706	Dist
C(5) C(6) C(4)	-3.049	0.620	Angle
C(18) O(19)	0.026	0.600	Dist
C(18) O(20)	-0.020	0.461	Dist

Angle at C(8) C(9) C(7)
Actual: 107.980 Hits: 20
Median: 111.597 Stdev: 3.585
Mean: 111.331
Min: 101.773 Max: 116.920
Lower Q: 110.289 Upper Q: 114.096
Z-score: -1.00

Examine in Mogul | Restrain
Restraining Everything

Cancel

Each bond length and angle is assigned a z -score:

$$z = \frac{\text{reported} - \text{CSDideal}}{\sigma(\text{CSDideal})}$$



Mogul Results

Results | Atoms | Z Distribution | Information

Z-scores: zero is best agreement. See information tab.
Click a column header to sort data.

Query	Delta	\sqrt{Z}	Type
C(8) C(9) C(7)	-3.617	1.009	Angle
C(16) C(17) C(15)	-3.419	0.810	Angle
C(4) C(5) C(3)	-2.681	0.717	Angle
C(4) C(5)	0.030	0.706	Dist
C(5) C(6) C(4)	-3.049	0.620	Angle
C(18) O(19)	0.026	0.600	Dist
C(18) O(20)	-0.020	0.461	Dist

Angle at C(8) C(9) C(7)

Actual:	107.980	Hits:	20
Median:	111.597	Stdev:	3.585
Mean:	111.331		
Min:	101.773	Max:	116.920
Lower Q:	110.289	Upper Q:	114.096
Z-score:	-1.00		

Examine in Mogul Restrain

Restraining Everything

Cancel

$$z = \frac{\text{reported} - \text{CSDideal}}{\sigma(\text{CSDideal})}$$




MOGUL and CRYSTALS

Simple interaction

- Finding out information about an odd looking bond and valence angle.

Richard Cooper
www.xtl.ox.ac.uk





MOGUL and CRYSTALS

Combined figures of merit.

A given atom is generally involved in several bonds and angles, each with an associated z-value.



$$\text{FOM} = \prod z_i$$

Individual z_i values are limited to the range 1 to 10.



MOGUL and CRYSTALS

Automatic interaction:



Traps incorrect structures.



What if the fragment does not exist in the CSD?

An interesting observation...



- Decide if the crystallography is secure.
- If it is, publish for the benefit of the community!

Validating the CSD

Mogul can be used to validate the structures stored in the CSD.


700 structures were extracted from the CSD and run through MOGUL.

The following are examples for which the FOMs were anomalous, possibly pointing to incorrect structures.






2 hydrogen atoms
have fallen off.

Could be refinement
problem, or
transcription error.



Very rare environment.
Despite many queries
with no results, the
FOM fares well – each
atom occurs in at least
one fragment with a hit.





Poor isopropyl groups.

Monoclinic $P2_1$, $Z'=2$,
and $\beta=90.46^\circ$

Pseudo 2_1 operator along
x at $y=0$, $z=0.5$ relates
molecules.

SG might be $P2_12_12_1$.






Poor alkyl chains.


Triclinic P1, $Z' = 4$.

Csp^3-Csp^3 distance
of 1.12\AA .

Z distribution tail
very long.



Poor isopropyl group.
C-C bonds 1.34 Å and
1.45 Å [ideal is 1.52 Å].
C-C-C angle is 104.3°
[ideal is 110.4°]



Z distribution shows
these outliers. Possible
outliers in F_0 data?

